



UnderDuct-VE Temperature Test Results:

<u>10-min. Test</u>	<u>200-F at 10-min.</u>	<u>300-F at 10 min.</u>	<u>400-F at 10 min.</u>	<u>750-F at 10-min.</u>
HEAT	no change	no change	color change only	color change only
ASTM #D2412	262 lbs.	320 lbs.	340 lbs.	270 lbs.
AS PSI	72.77	88.88	94.44	75

MONOXIVENT
Structural Capillary Systems
 PROPRIETARY AND CONFIDENTIAL
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MONOXIVENT. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF MONOXIVENT IS PROHIBITED.

	NAME	DATE	<i>Monoxivent</i> 1306 Mill St., Rock Island, IL 61201	REV
DRAWN	BP	10/15/2014	TITLE:	Duct Temperature Testing Setup
REV 1			FILE NAME	DTT Assembly
REV 2			SIZE: B	SCALE: 1:12
REV 3				SHEET 1 OF 1
REV 4				
REV 5				
REV 6				
REV 7				
REV 8				